

True Multi-Specialty
Portable Color
Veterinary
Ultrasound System



SONOVET R3



see it all ● ● ●
MEDISON

Ideally packaged for various animal veterinary applications, SONOVET R3 features color and power Doppler, Freehand 3D and cardiac imaging capabilities. And thanks to its compact design, it is designed to be with you wherever there is a power outlet. SONOVET R3 will provide you with optimal coverage in all fields of various animal veterinary.



Full-featured digital imaging on your desktop

The SONOVET R3 is a full-featured portable digital color ultrasound system with virtually all the imaging capabilities of a cart-based system.

- Full Spectrum Imaging™
- Harmonic imaging
- Color, power and PW Doppler
- Quick Scan™
- Freehand 3D
- Cardiac package

Organize and integrate your workflow

While the SONOVET R3 is a versatile diagnostic tool in its own right, it is an equally powerful collaborative tool. In addition to the basic review, analysis, networking, and archival capabilities of the SonoView™ provides the most advanced review and analysis capabilities available.

- SonoView™ Pro
- DICOM 3.0
- USB 2.0
- Probe: C2-4/20, CN2-8, CN4-9, L5-12/60, LN5-12/40, LV2-5/180, LV5-12/60

Expand your imaging horizons

Last but not least, the SONOVET R3 delivers powerful ultrasound imaging in an ergonomic package that can be used anywhere. The intuitive control interface features a full keyboard and large, backlit keys for all major functions. A custom cart featuring probe holders and plenty of mounting room for peripherals is also available.


- 15" LCD monitor
- Two probe ports
- Ergonomic handle in front panel
- Easily height adjustable control panel
- Lightweight



Abdomen vascular in 2D



Kidney in 2D



Mitral valve in pw Doppler