



The right probe for any application

The SONOACE X8 optional probes provide exceptional imaging performance in all applications. Advanced technology of automatic 3D volume reconstructing and expanded functional bandwidth delivers increased amount of clinical information and qualified diagnosis.

SONOACE X8

www.medison.com | info@medison.com

SONOACE X8



see it all
MEDISON

PG-MD-X8-CM-09.01.2008

SONOACE X8
























Probe Guide

Version 1.0

see it all
MEDISON

Probe Guide

																							
TRANSDUCER	C2-5EL	C2-5EP	C3-7EP	NEV4-9ES	NER4-9ES	EV4-9ED	ER4-9ED	L5-12EC	HL5-12ED	L5-12/50EP	L5-12/40EP	L4-7EL	C4-9/10EDN	3D2-6ET	3D4-8ET	3D4-8EK	3D5-9EK	P2-4AA	P3-7AC	P2-4AH	P3-5AC	2.0CW	4.0CW
Array Type	Curved	Curved	Curved	Curved	Curved	Curved	Curved	Linear	Linear	Linear	Linear	Linear	Curved	Curved	Curved	Curved	Curved	Phased	Phased	Phased	Phased	CW	CW
Sector Angle/ Field of View	76°	75.5°	70°	150°	150°	148°	148°	40mm	40mm	50mm	40mm	38mm	150°	84°	84°	69°	128°	16mm	9.6mm	19mm	14mm	N/A	N/A
Radius of Curvature	40mm	40mm	50mm	10mm	10mm	10mm	10mm	Flat	Flat	Flat	Flat	Flat	10mm	40mm	40mm	40mm	10mm	Flat	Flat	Flat	Flat	Flat	N/A
Broadband Frequency Range	2MHz~5MHz	2MHz~5MHz	3MHz~7MHz	4MHz~9MHz	4MHz~9MHz	4MHz~9MHz	4MHz~9MHz	5MHz~12MHz	5MHz~12MHz	5MHz~12MHz	5MHz~12MHz	4MHz~7MHz	4MHz~9MHz	2MHz~6MHz	4MHz~8MHz	4MHz~8MHz	5MHz~9MHz	2MHz~4MHz	3MHz~7MHz	2MHz~4MHz	3MHz~5MHz	2MHz	4MHz
APPLICATION																							
Abdomen	●	●	●											●	●	●							
Gynecology	●	●	●	●	●	●	●							●	●	●	●						
OB	●	●	●	●	●	●	●							●	●	●	●						
OB Early	●	●	●	●	●	●	●							●	●	●	●						
Fetal Heart	●	●	●	●	●	●	●							●	●	●							
Musculoskeletal								●	●	●	●	●											
Urology				●	●	●	●										●						
Pediatric								●	●	●	●	●	●					●	●	●			
Small Part								●	●	●	●	●	●									●	
Breast								●	●	●	●	●											
Vascular								●	●	●	●	●	●										
Cardiology																		●	●	●	●	●	●
TCD																		●	●	●	●		
Neonatal													●					●	●	●	●		
NEEDLE GUIDE	DISPOSABLE	N/A	REUSABLE	REUSABLE/ DISPOSABLE	REUSABLE/ DISPOSABLE	N/A	N/A	DISPOSABLE	DISPOSABLE	DISPOSABLE	N/A	DISPOSABLE	N/A	DISPOSABLE	DISPOSABLE	REUSABLE	REUSABLE	N/A	N/A	N/A	N/A	N/A	N/A
Biopsy Kit	BP-KIT-010		BP-C3-7IM	REUSABLE BP-KIT-019 DISPOSABLE BP-KIT-018	REUSABLE BP-KIT-019 DISPOSABLE BP-KIT-018			BPL-75/40AC	BP-KIT-022	BP-KIT-020		BP-KIT-014		BP-KIT-017	BP-KIT-017	VF-BP-PEC51	VF-BP-PEC47						
	